

ACCU-CHEK®
SmartGuide



**MOBILE DEVICE
AND SMARTWATCH
COMPATIBILITY**

繁體中文

行動裝置相容性

Accu-Chek SmartGuide app羅氏全智航應用程式和Accu-Chek SmartGuide Predict app羅氏全智航預測應用程式與滿足下列系統要求的大多數行動裝置相容。在您的行動裝置上安裝及操作應用程式前，請先檢查這些系統要求。

- 您的行動裝置必須支援Bluetooth藍牙低功耗技術5.0或更高版本。
- 您的行動裝置必須執行下列其中一個作業系統（OS）版本（不支援測試版作業系統）：
 - 對於iPhone，iOS版本必須是iOS 15.3或更高版本。
 - 對於Android裝置，作業系統版本必須是Android 11或更高版本。

智慧型手錶相容性

如果您使用的是iPhone，則可以將Accu-Chek SmartGuide app羅氏全智航應用程式與Apple Watch結合使用。

對於Apple Watch，watchOS版本必須是watchOS 9或更高版本。

其他智慧型手錶與應用程式不相容。

印刷版

如果您需要本文件的印刷版，請聯絡糖尿病照護部門/客戶服務部。印刷版免費提供，將在數天內寄送給您。

下載

在連線到網際網路時下載本文件並將其儲存到您的行動裝置，以便在沒有網際網路連線的情況下使用。本文件可從以下網址下載：www.accu-chek.com。*

* 下載可能會使用您的數據流量/產生費用。

ACCU-CHEK and ACCU-CHEK SMARTGUIDE are trademarks of Roche.
Apple Watch, watchOS, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries.
App Store is a service mark of Apple Inc., registered in the U.S. and other countries.
IOS is a trademark or registered trademark of Cisco in the U.S. and other countries.
Android, Google Play, and the Google Play logo are trademarks of Google LLC.
The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Roche is under license.
All other product names and trademarks are the property of their respective owners.

© 2025 Roche Diabetes Care



Roche Diabetes Care GmbH
Sandhofer Strasse 116
68305 Mannheim, Germany

www.accu-chek.com



Last update: 2024-12
1000090903(01)